



## RX400M

### The Medical USB Solution - fast, easy & secure

Save your entire collection of diagnostic data in full resolution with just the touch of a finger. Deliver studies on a state of the art write protected Patient USB-Stick with the RX400M high-performance data copy system. The USB Sticks can be used alongside conventional patient CD, with added benefits such as high speed, storage space and independence from the availability of a CD drive.

The compact RX400M is easy to set up for system integrators, health care providers and hospitals.

#### Functionality

The RX400M is a compact, fully integrated desktop device including integrated computer, 4 front USB ports, label printer and touch screen and can be complemented with an internal DVD drive to import data from patient CD/DVD into PACS. Combined with leading DICOM solutions the system can be easily integrated into your existing workflows like a DICOM printer. Patient sticks are produced as easy as you know it from patient CDs with the advantages such as fast opening of data and high storage capacity that allows multiple studies.

The Rimage Patient USB-Sticks are automatically write-protected so it can be used safely, and computer viruses don't stand a chance. The integrated printer prints the patient label as soon as the finished stick is unplugged.

#### What makes a Patient USB-Stick successful?

With the Rimage patient stick we enhanced data security with unalterable data. This is achieved by WORM procedure. With "Write Once, Read Many" data on the Rimage patient USB-Stick cannot be modified. This provides high security against virus and malware, as the sticks can only be written by the RX400M system. Other USB-Ports can only read but not write any data on Rimage patient USB-Sticks.

With our Top Line and BaseLine Pro Stick, the Rimage RX400M system will be able to add more data to your existing data and protect the USB-Stick afterwards with a new write protection. Here are some further details of the Rimage Patient USB-Stick:

#### OVERVIEW RIMAGE PATIENT USB-STICKS

##### TopLine Patient USB-Stick (64GB)

##### BaseLine Pro Patient USB-Stick (16GB)

USB Version: 3.1

Multi-Write Protect Feature

Connector: USB A

Max Read/Write Speed: 400/200 MB/s

Sequential read/write: 40/30 MB/s

Data Retention Time: up to 10 years

##### BaseLine Patient USB-Stick (16GB)

USB Version: 2.0

SingleWORM Feature

Connector: USB A

Max Read/Write Speed: 200/100 MB/s

Sequential read/write: 25/8 MB/s

Data Retention Time: up to 3 years



## Value for Hospitals and Clinics

---

### Simplified Operations

The integrated touchscreen shows all waiting studies, sorted by patient, coming from any connected DICOM source in your hospital. Via touch screen, you can select which study you want to store on the USB-Stick.

With the high-performance hardware you can create USB-Sticks with common data, approx. 500 MB, in less than 30 sec. This allows you to create media on the spot, whenever a patient shows up at your front desk.

### Security

The RX400M can exclusively be used with Rimage USB-Sticks that have special security features such as WORM (write once, read many) protection. This dramatically increase the propability that no other files that may contain malware, viruses or other threats can be added on the stick and thus get into highly sensitive networks. Only with the RX400M system you are able to unlock the pre-protection, add new data, and activate the WORM protection or multi-write protection mode afterwards. This limits the risks of Malicious USB-Sticks and provides an additional IT-Environment protection.

With the integrated, adhesive paper thermal printer you can easily label the USB-Sticks individually, ensuring you are handing the right stick to the right person. The dedicated patient labels are automatically created on a PE-Foil label once you unplug a recorded and verified USB-Stick.

### Reliable Workflow

Without any moving and mechanical parts, the RX400M is very service friendly. Most of the parts are easily swappable and accessible by opening the cover.

With the easily accessible SSD that is securely locked within an cage, containing device configuration data and locally used study data, you can even request an exchange unit ifquired and just use the existing SSD in the new unit. In this case, your sensitive data will never leave your environment.

## CONTACT INFORMATION



**Rimage Corporate Headquarters**  
**201 General Mills Blvd**  
**Golden Valley, MN 55426**  
**USA**



**+1.800.445.8288**



**Sales@rimage.com**

**R I M A G E® RX400M**

## Value for System Integrators

---

### Dedicated Support

The Rimage worldwide service organization has dedicated applications engineers to provide support during all stages of system integration.

### Easy Integration

The RX400M API enables easy integration for DICOM Software providers. The RX400M is already fully integrated into several DICOM patient CD/USB Solutions.

### USB Technology

Allows Integrators to easily output on a new media source besides the Rimage optical systems.

## Value for Patients

---

### Reliable Records

Having the option of reusing your USB-Stick, you can keep track of all your medical data on a single source. This also allows doctors to see the entire patient's history, which helps them to better diagnose disease progression.

### Shorter Wait Times

RX400M will provide a USB-Stick with a typical study in about 30 seconds, so patients don't have to wait long for their personalized USB-Sticks.

### Data Privacy and Security

With the RX400M we are adding a new layer of security for patient records on physical data. Your personal medical data can be encrypted and in case of a loss, it will be nearly impossible to access your sensitive patient data and records. This feature however, needs to be supported by your DICOM Software.

## Service Options

---

For the RX400M system we offer two different types of maintenance contracts to accommodate the needs of our customers.

### Send-in Service

Customer needs to send in the RX400M system at his own costs. When the system arrives in the service depot, the Service team will repair the system as quickly as possible. After repair, the system will be sent back to the customer and Rimage is covering the return freight costs.

### Fast-Track Service

Rimage will send out a replacement unit to the customer on next working day on Rimage costs. When the system arrives at the customer, the customer needs to send back the broken unit using the included return label.